



## **BipoTrode**

### 5 Fr. Bipolar Electrode

5 Fr. bipolar dissection, preparation and vaporisation • Can be used with common hysteroscopes with a 5 Fr. working channel • Autoclavable • Robust design

# BipoTrode

## 5 Fr. Bipolar Electrode



- For precise dissection, preparation and coagulation using saline distention solutions
- Small diameter of 5 Fr. allows use with hysteroscopes
- Made of semiflexible stainless steel it is suitable for use with curved working channels and ideal for use through the straight working channels of the compact R. WOLF hysteroscopes
- Robust and ergonomic handle
- Autoclavable and reusable

**BipoTrode bipolar electrode,**

5 Fr., WL 360 mm .....823005

also:

**HF bipolar connecting cables, 3 m**

for

WOLF / Martin / Berchtold / Aesculap .....8108.131

ERBE / ACC / ICC / T .....8108.132

2 pin Valleylab .....8108.133

Recommended for the use with the following R. WOLF hysteroscopes:

**5 mm continuous-irrigation outpatient operating hysteroscope, 20°**

with continuous-irrigation, by Grochmal,

working channel 5 Fr., WL 217 mm .....8986.401

**4 mm mini compact continuous-flow operating hysteroscope, 12°**

with oblique eyepiece and straight working channel with

continuous irrigation, working channel 5 Fr.,

WL 230 mm .....8977.411

**5,5 mm *OptiFLOW*™ continuous-flow operating sheath,**

working channel 5 Fr. ....8974.351

works with:

**PANOVIEW telescope**

∅ 2,7 mm, WL 300 mm

direction of view 30° .....8974.402

12° .....8974.412